


ASPIRATORI AD ANELLO LIQUIDO NON CARENATI

Accessori

	Kw	Volt	Hz	g/min	Vuoto: mbar	Portata lt/min	Dim		dB(A)	ITALIANO	ENGLISH	FRANCAIS	ESPAÑOL
P001 * a richiesta anche a 60Hz e trifase	1,5	230	50	2800	558,6	750	40x32x42 40Kg	4	66	-centralina -filtro secreti -elettrovalvola -manometro -rid. pressione -valv. reg. vuoto -vuotometro	-control board -secretions filter -electrovalve -pressure gauge -pressure reducer -vacuum reg. valve -vacuum meter	-Boîte de contrôle -Filtre sécrétions -électrovanne -manomètre -réducteur de pression -clapet de contrôle du vide -vacuomètre.	-Con Centralita -Filtro de secreciones -Electroválvula -Manometro -Reductor de presión -Válvula reg. de vacío -Vacúmetro
P001/S	1,5	230	50	2800	558,6	750	40x32x42 40Kg	4	66	-come per P001 -vaso separatore	-like us P001 -air/water separator	-comme P001 -air/eau séparateur	-Como P001 -Vaso separador
P001/SB	1,5	230	50	2800	558,6	750	55x42x46 45Kg	4	66	-come P001 -vaso separatore -base metallica	-like us P001 -air/water separator -metal support	-comme P001 -air/eau séparateur -support en métal	-Como P001 -Vaso separador -Base metálica
P001/P	1,5	230	50	2800	558,6	750	40x32x42 40Kg	4	63	-come P001 -polmone silenziatore	-like P001 -silencer	-comme P001 -silencieux	-Como P001 -Silenciador
P001/T	1,5	230	50	2800	558,6	750	40x32x42 40Kg	4	63	-come P001 -vaso separatore -polmone silenziatore	-like us P001 -air/water separator -silencer	-comme P001 -air/eau séparateur -support en métal	-Como P001 -Vaso separador -Silenciador
P001/TB	1,5	230	50	2800	558,6	750	55x42x44 45Kg	4	63	-come P001 -vaso separatore -polmone silenziatore	-like us P001 -air/water separator -silencer -metal support	-comme P001 -air/eau séparateur -silencieux -support en métal	-Como P001 -Vaso separador -Silenciador -Base metálica
P002 * a richiesta anche a 60Hz	0,75	230	50	2800	532	383	32x27x37 35Kg	2	59	-centralina -filtro secreti -elettrovalvola -manometro -rid. pressione -valvola reg. vuoto -vuotometro	-control board -secretions filter -electrovalve -pressure gauge -pressure reducer -vacuum reg. valve -vacuum meter	-Boîte de contrôle -Filtre sécrétions -électrovanne -manomètre -réducteur de pression -clapet de contrôle du vide -vacuomètre.	-Con Centralita -Filtro de secreciones -Electroválvula -Manometro -Reductor de presión -Válvula reg. de vacío -Vacúmetro
P002/S	0,75	230	50	2800	532	383	32x27x37 35Kg	2	59	-come P002 -vaso separatore	-like us P002 -air/water separator	-comme P002 -air/eau séparateur	-Como P002 -Vaso separado



P001 - P001/S - P001/P - P001/T




P001/SB - P001/PB - P001/TP



P002 - P002/S

ASPIRATORI AD ANELLO LIQUIDO CARENATI

Accessori

	Kw	Volt	Hz	g/min	Vuoto: mbar	Portata lt/min	Dim		dB(A)	ITALIANO	ENGLISH	FRANCAIS	ESPAÑOL
UCA1 * a richiesta anche a 60Hz e trifase * a richiesta kit super insonorizzante	1,5	230	50	2800	558,6	750	55x42x46 57Kg	4	61	-come P001 -vaso separatore -carena metal	-like us P001 -air/water separator -metal cabinet	-comme P001 -air/eau séparateur -support en métal	-Como P001 -Vaso separador -Carenado metálico
UCA2	0,75	230	50	2800	532	383	40x32x46 45Kg	2	56	-come P002 -vaso separatore -carenatura metal	-like us P002 -air/water separator -metal cabinet	-comme P002 -air/eau séparateur -support en métal	-Como P002 -Vaso separador -Carenado metálico



UCA1



UCA2