

SEALER SL13

NEW



GUARANTEED TO MAINTAIN STERILIZATION PERFORMANCE

SL13 an **electronically controlled thermo-sealer** with a sealing band of 12 mm fit to resist the sterilization processes of the modern sterilizers featuring **fractioned pre-vacuum**, and to assure the sterility of the packages.

CONSTANT HIGH PERFORMANCE

Efficient heat control in the welding area is achieved thanks to the **automatic temperature adjustment**, ensuring high performances. **SL13** also feature an overheating protection device and an automatic **energy-saving standby system**.

ERGONOMIC AND PRACTICAL DESIGN

The **solid structure** in plastic and aluminium and the **ergonomic and minimalist shape** help to maintain a high level of hygiene, weight and overall size and **excellent stability** during use. All this with the objective of offering a product suited to the daily needs in the medical area.

TECHNICAL FEATURES

GENERAL FEATURES

Power supply voltage	220/240 Vac - 50/60 Hz
Nominal Power	150 W - 0.45 A
External dimensions (without roll)	474 x 374 x 200 mm (LxDxH)
Net Weight	5 kg
Type of sealing	Continuous cycle and constant temperature
Nominal operating temperature of the heating resistance	180 °C
Max sealing length	12 mm
Maximum welding width	300 mm
Sealable material (in accordance to the Directive EN 868)	Paper /polypropylene film
Max diameter of the roll	200 mm



Application of the Directives: 2006/95/EC - 2004/108/EC - Reference Standards: EN 61010-1:2001, EN 61326-1:2006. 24-months guarantee.



SINCE 1948: EXPERIENCE
AND INNOVATION

FARO S.p.A.

via Faro, 15 - 20876 Ornago (MB) - Italy Tel. +39 039.68781 - Fax +39 039.6010540
www.faro.it - comm.estero@faro.it

FARO FRANCE

Za Tgy Coriolis - 71210 Monchanin - France Tel. +33 385.779680 - Fax +33 385.779688
www.farofrance.com - farofrance@farofrance.com

FARO DEUTSCHLAND GMBH

Gewerbepark Heideckhof Heideckstr. 179 D-47805 Krefeld - Germany Tel. +49 2151.936921 - Fax +49 2151.936933
www.farodeutschland.de - info@farodeutschland.de

FARO IBERICA Medicodental S.L.

C/ Tarragona, 16 local - 43830 - Torredembarra (Tarragona) - Spain Tel. +34 9776.41901 - Fax +34 9771.3085
www.faro.it - info@faroiberica.com

Certified Company



ISO 9001:2008 ISO 13485:2003

INFO LINE:

